

Pay As You Go: Self Funded Schedule of Fees

Service Category	Items of Service	Prices	Totals per Category
Free Hearing Check (Screen)	Otoscopy 4 Frequency Screen @ 250Hz, 500Hz, 1KHz & 2KHz – (determines whether further testing/referral required)	\$0	
Full Diagnostic Hearing Test	Medical & Audiological History Taking & Communication (COSI) Goal setting Documentation	\$20	
	Otoscopy/Video Otoscopy	\$20	
	Tympanometry (Bilateral)	\$20	
	Wide Band Pure Tone Audiometry 250Hz- 10KHz	\$40	
	Speech Discrimination Noise Test Battery: SPIN, AB Words, Quick Sin	\$20	\$120
<i>To Acquire</i>	<i>Oto Acoustic Emissions (AOE's) Test for Outer Hair Cell Function</i>	<i>\$20</i>	
Report to G.P.	Audiological Report (Take Home) & Letter Send to G.P./Specialist (if required)		\$60
	Hearing Aid recommendation/Counselling (As Required) - full hearing aid demonstration		
			Total \$180
Hearing Aid Fittings And Adjustments	Hearing Aid Programming & Fine Tuning (based on client feedback) - Installation of additional programs - Connection & Calibration of accessories e.g. remote control, Streamer, remote microphone etc. - Aided Speech discrimination in noise test (e.g. SPIN , Quick Sin or Live Voice Sentences) Verification of Hearing Aid Fitting (Probe Tube REM's) correctly calibrated within 5mm of T.M. (Need to Include to show professionalism and verify aid output)	80 60	 \$140
	Connection with Smart Phone App or accessory and Coaching in Use of App/multiple programs etc.	20	\$160
1 Follow Up Minimum Recommendation	30 minutes. To review of COSI Goals, Hearing Aid Performance and acclimatisation, Physical and Acoustic comfort and fine tuning adjustment (as required)	\$60	
2nd Three Month Adaptation Review	30 mins. Check usage, Review acclimatisation.	\$60	
Specialised Ear Wax Removal	- Ear Wash Lavage System - Addition of Melaleuca oil (Anti Bacterial/Anti-flamatory) And/OR	\$50 per ear	\$100
	Unique Dry Combination Technique - Micro Suctioning (using Microscope) - Use of Curettes & Forceps (As required under microscope)	\$50 per ear	\$100